

Electrical Stimulation Machine Comparison Table

	NEUROTRAC PELVITONE	NEUROTRAC CONTINENCE	NEUROTRAC TENS	ACCURATE TENS N602	NEUROTRAC OBSTETRIC TENS	NEUROTRAC SPORTS
TENS	IN CUSTOM	IN CUSTOM	YES	YES	YES	Prog 1, 2, 12 only
NMES	YES	YES	NO	NO	NO	YES
IFT	NO	NO	NO	NO	NO	NO
CUSTOM PROGRAMS	3	4	3	YES – ALL	NO	3
PRESET PROGRAMS	11	9	11	NO	YES:ONE	15
mA Range	0 - 90	0 - 90	0 - 80	0 – 90 10mA/unit	0 - 80	0 - 90
Hz Range	2 - 100	2 - 100	2 - 200	2 - 150	150 bursts rest / pre-set 90 boost	2 - 100
uS Range	50 - 330	50 - 450	50 - 300	50 - 300	200 Cont mode, 220 Boost mode	50 - 450
Minutes Range	1 - 90 min	1 - 90 min	1 min - 12 hours	30 min, 60 min, Continuous	Continuous	1 - 90 min
Work Time Range	2 - 90 sec	2 - 90 sec	N/A	N/A	N/A	2 - 99 sec
Rest Time Range	2 - 90 sec	2 - 90 sec	N/A	N/A	N/A	2 - 99 sec
Ramp Up	0.1 - 9.9 sec	0.1 - 9.9 sec	N/A	N/A	0.8 Rest to Boost	0.1 - 9.9 sec
Ramp Down	NO	0.1 - 9.9 sec	N/A	N/A	N/A	NO
Alter SYN/ALT	YES	YES	NO	NO	N/A	YES
Set Delay Between Channels	YES	YES	NO	NO	N/A	YES
LOCK Function	YES	YES	YES	NO	NO	YES
Average mA used display	YES	YES	YES	NO	NO	YES
Total hours used display	YES	YES	YES	NO	NO	YES
Continuous	Custom prog	Custom prog	YES	YES	YES	Not in custom
Modulate function	NO	NO	YES High/Low	YES uS or uS+Hz	YES in preset	Only some pre-sets